



# Longitudinal Imaging Biomarkers of Disease Progression in DLB

---

BIOSPECIMEN COLLECTION & PROCESSING

# Overview

---

1. Specimen uniformity and quality
2. Site Equipment
3. Procedures
  - Kit Contents and Ordering
  - Sample Labelling
  - Sample Collection and Processing
  - Shipping Samples
  - Non-Conformance
4. Contact Information

# Specimen Uniformity and Quality

---

GENERAL REMINDERS

# Specimen Standardization and Quality

---

Most biomarkers are sensitive to *time* and *temperature*

- Standardization of processing across sites is key
- Specimens must be processed within 2 hours of collection
- Reference the *BioSEND Biomarker Specimen Collection, Processing, and Shipment Manual* as needed
- Do not replace or supplement any kit components without first receiving approval from BioSEND/NINDS

Questions? Email [biosend@iu.edu](mailto:biosend@iu.edu)

# Site Consumables and Equipment

---

*Sites will need to supply the following items:*

- Gloves
- Alcohol wipes
- Butterfly needles
- Tourniquet
- Gauze pads
- Bandages
- Microcentrifuge tube rack
- Sharps bin and lid
- Crushed ice
- Pipettes and pipette tips
- Centrifuge capable of maintaining 4°C
- -80°C Freezer
- Dry ice

# Procedures

---

MAINTAINING SPECIMEN UNIFORMITY AND QUALITY

# Biospecimen Collection Protocol

[illegible]

# Kit Contents and Ordering

---

- All sites will be sent a Supplemental Kit with their first kit shipment
  - Contains extra blood collection tubes and processing supplies
  - May be used to replace items in study visit kits
- Study Visit Kits should be ordered as soon as visits are planned
  - Contains collection, processing, and shipping supplies specific to each visit
  - Includes barcoded labels
  - The supplies/labels in each study visit kit are intended for that visit only




# Kit Contents and Ordering – REDCap Survey

<http://kits.iu.edu/biosend/dlb>

Order kits online through the Kit Request Module for:

- Blood kits
- Supplemental Kit
- Extra Supplies

Please provide as much notice as possible when ordering kits and/or supplies.



The screenshot shows the BioSEND NINDS Biosend DLB Kit Request System interface. At the top, the logo "BioSEND" is displayed in large blue letters, with "Biospecimen Exchange for Neurological Disorders" written below it. The title "NINDS Biosend DLB Kit Request System" is centered. Below the title, there is a "Study Site" label and a dropdown menu. A red asterisk and the text "\* must provide value" are positioned below the "Study Site" label. A "Submit" button is located at the bottom right of the form.

**BioSEND**  
Biospecimen Exchange for Neurological Disorders  
**NINDS Biosend DLB Kit Request System**

Study Site   
\* must provide value

**Submit**

# Kit Contents and Ordering: Confirm Site Info

DLB Kit Request Module

<b>Study Site</b> <small>* must provide value</small>	<input type="text" value="Indiana University"/>	<div>Select your site from the drop-down list</div>
<b>Indiana University School of Medicine</b>  Carolyn Dunifon Dept. of Medical & Molecular Genetics 351 West 10th Street, TK-318 Indianapolis, IN 46202-3002 (317) 274-5751 cdunifon@iu.edu		
<b>Is the contact name above correct?</b> <small>* must provide value</small>		<div><input type="radio"/> Yes <input checked="" type="radio"/> No</div> <div>reset</div>
<b>New Contact Name</b> <small>* must provide value</small>		<input type="text" value="Claire Wegel"/>
<b>Is the shipping address above correct?</b> <small>* must provide value</small>		<div><input type="radio"/> Yes <input type="radio"/> No</div> <div>reset</div>
<b>Is the e-mail address above correct?</b> <small>* must provide value</small>		<div><input type="radio"/> Yes <input type="radio"/> No</div> <div>reset</div>
<b>Is the phone number above correct?</b> <small>* must provide value</small>		<div><input type="radio"/> Yes <input type="radio"/> No</div> <div>reset</div>

# Kit Contents and Ordering: Kit Types

## DLB Kit Request Module

- Kits are not specific to a subject or time point. After collection, sites will indicate the subject and time point to which BioSEND should link the samples.

<b>Kit Type</b> **Please allow two weeks for shipment** * must provide value	<div><input checked="" type="checkbox"/> Baseline or Annual Visit Kit</div> <div><input type="checkbox"/> Supplemental Kit</div> <div><input type="checkbox"/> Extra Supplies</div> <p>Please specify in comments if you need kits before the standard two week shipment time.</p>
<b>CSF Sprotte® Needle Gauge</b> * must provide value	<div>22</div> <div>24</div> <div>reset</div>
<b>Baseline or Annual Visit Kit Quantity</b> * must provide value	<div>2</div>

# Kit Contents and Ordering: Kit Breakdown

DLB Kit Request Module

Comments

Expand

Each Baseline or Annual Collection Kit

Blood Collection Kit Contents:

2 - Lavender-top EDTA tube (10 ml), glass  
2 - Purple-top EDTA tube (3 ml), plastic  
2 - Red-top serum tube (10 ml), glass  
2 - PAXgene® tube (2.5 ml)  
15 - Siliconized cryovials (2ml)  
2 - Disposable transfer pipette (3ml)  
2 - 15ml conical tube (individually wrapped)

CSF Collection Kit Contents:

11 - Siliconized cryovials (2ml)  
2 - 15ml conical tube (individually wrapped)  
2 - 50ml conical tube (individually wrapped)  
1 - Medication Transfer Filter Straw  
1 - LP tray with 22 gauge Sprotte® needle

1 - 25-slot cryobox  
1 - Shipping label packet  
1 - Airway bill envelope  
1 - Shipping container for dry ice shipments  
2 - Biohazard bag with absorbent sheet  
8 - individual tube bubble pouch  
30 - Cryohold pre-printed labels

Submit

Kit contents of selected kit will appear at the bottom of the page


# Automatic Kit Shipments

---

- After subject completes baseline visit and BioSEND receives BL samples, BioSEND sets up automated kit sending schedule for subject's subsequent visits
- Schedule gives 2 month window around the longitudinal study visit target (1 month on either side)
- BioSEND will send kits prior to start of study window
  - Reduces effort for study coordinators
  - Sites only need to order kits if visit will occur AHEAD of the study visit window
- All study visit target dates are determined from Baseline Visit (not from last study visit date)

# Automatic Kit Shipments

---



# Kit Contents and Ordering: Blood Kit


---





# Kit Contents and Ordering: CSF Supplies

CSF:



LP Tray:





# Collection Volumes

Total blood and CSF volumes

Sample Type	Amount
Whole Blood for RNA	5 ml
Whole Blood for Plasma and Buffy Coat	20 ml
Whole Blood for Serum	20 ml
Whole Blood for Banking	6 ml
Cerebrospinal Fluid	10 ml

# Kit and Supply Ordering

---


- Click “Submit” to send order to BioSEND; staff will confirm receipt of your order
- Please allow two week turn-around time for kit shipments
- If urgent request needed, please note date needed by in comments and email BioSEND. We cannot guarantee urgent orders, but we will do our best to accommodate.
- BioSEND will send confirmation of shipment and tracking number when supplies ship

# Sample Labelling: Example Labels

---

## Labels are provided by Indiana University

- Please check that all samples are properly labelled to ensure correct identification by IU
- If do not have enough labels to complete a visit, please contact IU *immediately*
- Labelling the tubes during processing prevents sample mix-ups




# Sample Labelling: Label Placement

---


Please...

- Label all collection and aliquot tubes before cooling, collecting, processing, or freezing samples
- Label only 1 subject's tubes at a time to avoid mix-ups
- Wrap the label around the tube horizontally - label position is important for all tube types
- Make sure the label is completely adhered by rolling between your fingers




# Sample Collection and Processing


Blood Tube Draw Order




# Sample Collection and Processing: Serum



# Sample Collection and Processing: Whole blood RNA




# Sample Collection and Processing: Plasma & Buffy Coat






# Sample Collection and Processing: Whole Blood

---



# Sample Collection and Processing: CSF




# Sample Collection and Processing: Aliquots

---


Filling biomarker serum, plasma, and CSF aliquots:

- Fill as many cryovials as possible to 1.0 ml (plasma, serum & CSF)
- Over-filled vials may burst in freezer!
- Ship ALL material to IU, even if final vial is less than standard volume



# Sample Collection and Processing: Timeline

Timeline for blood processing



# Sample Collection and Processing: Issue #1

---

## Troubleshooting Blood Collection

### Issue #1: Tube with little/no vacuum

- Always check expiration date on the tube before beginning blood draw and discard expired tubes
- Store tubes at “room temperature” – extreme temperature can affect vacuum
- Keep extra vacutainer tubes from supplemental kit nearby during blood draw to replace “bad” tubes
- If this is a frequent occurrence, report tube type and lot number to IU.

# Sample Collection and Processing: Issue #2

---

## Troubleshooting Blood Collection

### Issue #2: Hemolyzed serum and/or plasma caused by incorrect collection

Cause: Blood Collection Methods	Corrective Action
Improper venipuncture site	Draw from median cubital, basilic, and cephalic veins from antecubital region of arm
Prolonged tourniquet use	Tourniquet should be released after no more than 1 min, excessive fist clenching should be avoided
Not allowing alcohol to dry on skin before venipuncture	Without touching, allow the venipuncture site to air dry
Use of too large/small bore needle resulting in excess force applied to blood	Avoid using too small/large needle. Needle size dependent on the subject's physical characteristics & amount of blood to be drawn. Most commonly used sizes are 19 – 23.
Pulling/pushing plunger too fast while drawing/transferring blood	Avoid drawing the syringe plunger too forcefully when collecting blood
Ensure all blood collection assemblies are fitted securely, to avoid frothing	

For more information, visit: [http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk\\_Jan2004\\_VS7167.pdf](http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk_Jan2004_VS7167.pdf)

# Sample Collection and Processing: Issue #2 continued

---

## Troubleshooting Blood Collection

### Issue #2: Hemolyzed serum and/or plasma caused by incorrect processing

Cause: Sample Processing Methods	Corrective Actions
Vigorous mixing/shaking	Gently invert blood collection tube when mixing additive with specimen, follow guidelines in Biologics Manual regarding number of times to invert each type of tube
Not allowing serum to clot for recommended time	Serum tubes without clot activator should be allowed to clot for 60 min in a vertical position
Exposure to excessive heat or cold	Keep samples at ambient temperature until processing
Prolonged contact of serum/plasma with cells	Do not store uncentrifuged samples beyond recommended time

For more information, visit: [http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk\\_Jan2004\\_VS7167.pdf](http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk_Jan2004_VS7167.pdf)

# Sample Collection and Processing Form

Direct link:

<https://redcap.link/DLBSampleForm>

First part captures basic subject and visit information

# BioSEND

## Biospecimen Exchange for Neurological Disorders

Please verify/update the information below. When you click the "Submit" button below, a PDF copy of the Sample Record and Shipment Notification Form will be emailed to you.

Please print a copy of that document and **include it in the shipping container.**

### Longitudinal Imaging Biomarkers of Disease Progression in DLB Study

Study Site

Mayo Rochester ▼

Email address of staff member completing this form

*Note: A copy of the completed sample form and the shipping manifest will be sent to this address.*

GUID

Sex (used for DNA quality control)

Visit

Kit Number

Date of venipuncture blood collection

  Today M-D-Y

Time of venipuncture blood collection

  Now H:M

Use 24 Hour clock

 Returning?

AAA





# Sample Collection and Processing Form

Direct link:

<http://kits.iu.edu/biosend/DUBSampleForm>

Second part captures processing information

SERUM	
Number of SERUM aliquots shipped:	<input type="text"/> Each aliquot should be 1 mL
RNA PAXGENE	
Number of PAXGene™ tubes shipped:	<input type="text"/>
WHOLE BLOOD EDTA	
Number of WHOLE BLOOD tubes shipped:	<input type="text"/>
PLASMA EDTA	
Number of PLASMA EDTA aliquots shipped:	<input type="text"/> Each aliquot should be 1 mL
Number of BUFFY COAT aliquots shipped:	<input type="text"/>
CSF	
Was CSF collected at this visit?	<div><input type="button" value="Yes"/></div> <div><input type="button" value="No"/></div> <div>reset</div>
NOTES	
Please record any issues with collection/processing:	<div><input type="text"/></div> <div>Expand</div>

# Sample Collection and Processing Form

Direct link:

<http://kits.iu.edu/biosend/DLBSampleForm>

PDF form of responses will be emailed to you. **Print a copy of the Frozen Shipping Manifest and include with shipment.**

## DLB Mayo Sample Record and Shipment Notification Form <sup>Page 1</sup>

Please verify/update the information below. When you click the "Submit" button below, a PDF copy of the Sample Record and Shipment Notification Form will be emailed to you.

Please print a copy of that document and include it in the shipping container.

---

Longitudinal Imaging Biomarkers of Disease Progression in DLB Study

---

Study Site ☐ Mayo Rochester  
☐ Mayo Jacksonville

---

Email address of staff member completing this form

---

Note: A copy of the completed sample form and the shipping manifest will be sent to this address.

---

GUID

---

Sex (used for DNA quality control) ☐ Male  
☐ Female  
☐ Other

---

Visit ☐ BL  
☐ 12M  
☐ 24M  
☐ 36M  
☐ 48M  
☐ 60M  
☐ 72M  
☐ 84M  
☐ 96M  
☐ 108M

---

Kit Number

---

Date of venipuncture blood collection

---

Time of venipuncture blood collection

---

(Use 24 Hour clock)

---

# Sample Collection and Processing Form

Submission of the Shipping Manifest portion of the form serves as shipment notification to BioSEND

- Must be completed prior to shipment
- If samples are shipped and not received, BioSEND will follow-up with courier. It is recommended that sites also track shipment to ensure safe delivery

# Shipping Frozen Samples: Tips

---

## Packing and Shipping Frozen Samples

- Serum, plasma, buffy coats, CSF, whole blood and RNA all ship frozen
- Ship frozen samples on dry ice
- Frozen samples should be shipped **only** Monday through Wednesday
- Always fill carton to **top** with dry ice
- Do not pack shipment until the day of pickup



# Shipping Samples

## Packing and Shipping Frozen Samples

Class 9 Dry Ice Label should not be covered with other stickers and must be completed, or UPS will reject/return your package!

The image shows a Class 9 Dry Ice Label template. It features a large triangle with vertical bars inside, and a large '9' at the bottom. The label includes fields for shipper and consignee information, and a section for dry ice weight. Red arrows point to specific areas with labels: 'Net weight of dry ice in kg' points to the weight field, 'Your name & address' points to the shipper's name and address field, and 'IU information and address' points to the consignee's name and address field.

Shipper's Declaration not Required.

Dry Ice amount must be in kilograms.

Note: 2 lbs. = 1 kg.

Airwaybills / airbills must have the following:

1. Dry Ice; 9; UN 1845
2.  $\frac{\text{Number}}{\text{(Number pks)}} \times \frac{\text{wt}}{\text{(wt)}} \text{ Kg}$

Net weight of dry ice in kg

Dry Ice  kg.

Your name & address

Shipper's Name and Address

UN 1845

Consignee Name and Address

IU information and address

9

06426 1/01 RRD

# Shipping Samples: Frozen

Do not ship more than 4 biohazard bags in a single shipper (equivalent to two subject-visits).

- Allows room for dry ice to keep samples frozen in transit
- Minimizes loss in the rare but unfortunate event of courier issues



Shipping Samples – UPS: <https://kits.iu.edu/UPS>

[illegible]

# Shipping Samples via UPS

IU UPS ShipExec Shipping Portal

- Print out UPS air waybill
- Ensure all elements (barcode, return address, etc.) printed clearly
- Fold and insert UPS air waybill into clear plastic sleeve on package

JOHN SMITH INDIANA UNIVERSITY 410 WEST 10TH STREET INDIANAPOLIS IN 46202	2 LBS	1 OF 1
RS		
SHIP TO: SCHOOL OF MEDICINE 317-278-2694 INDIANA UNIVERSITY TK 217 351 W 10TH ST INDIANAPOLIS IN 46202		
	IN 461 9-01	
		
UPS NEXT DAY AIR		1
TRACKING #: 1Z 976 R8W 84 3985 8595		
 SAMPLE		
BILLING: P/P DESC: Biological Specimens RETURN SERVICE		
Reference No.1: 4087277		
XOL 20.03.09 NV45 83.0A 12/2019		 TM



# Shipping Samples: Closures

Date	Holiday
January 1	New Year's Day
3 <sup>rd</sup> Monday in January	Martin Luther King, Jr Day
4 <sup>th</sup> Monday in May	Memorial Day
June 19	Juneteenth (observed)
July 4	Independence Day (observed)
1 <sup>st</sup> Monday in September	Labor Day
4 <sup>th</sup> Thursday in November	Thanksgiving
4 <sup>th</sup> Friday in November	Friday after Thanksgiving
December 25	Christmas

Please also consider weather when shipping. UPS will post service updates on their webpage. Reach out to BioSEND if you are unsure if it is safe to ship.

# Non-Conformance

---

Non-conformance to standard procedures may reduce the utility of the biospecimens:

- Not processing serum/plasma within 2 hours of collection allows for breakdown of certain proteins and small molecules
- Over/under centrifuging changes plasma, serum, CSF composition



# Non-Conformance Reporting con't

---

Most common non-conformance issues:

- Shipment notification not sent
- Samples shipped for weekend/holiday delivery
- Sample form incomplete/inaccurate
- Low volume
- Unlabeled or mislabeled tube(s)
- Sample hemolysis



# Non-Conformance and Inventory Reporting

---

Most common non-conformance issues:

- BioSEND will notify sites directly of any issues upon receipt
- BioSEND will email sites a monthly inventory report of all samples received from that site to date
- ***If you are experiencing issues, please reach out to us for help! It is much easier to prevent an issue before sample collection & shipment than trying to fix it after the fact***

# BioSEND.org

On the website, you can:

- Access your study's kit request module and sample submission form
- Download the most recent version of the Manual of Procedures
- View a recording of this training
- Find information about holiday closures
- Access shipping resources

## Study Resources

### KIT REQUEST MODULE

Please follow the below link to access the Kit Request Module. This link will direct you to a REDCap database where study coordinators and staff may request kits, individual supplies, and/or labels. Please allow a total of two weeks for kit requests to be fulfilled.

[Kit Request System →](#)

### SPECIMEN COLLECTION AND PROCESSING FORM

Please use the below link to access the collection and processing form for this protocol. This form must be completed prior to shipment. We also ask that all shipments include a physical copy of the shipping manifest portion of the form.

[Specimen Collection and Processing Form →](#)

### MANUAL OF PROCEDURES

The below downloadable manual was created specifically for the DxCTEII study. Please feel free to explore the manual through the hyperlinked 'Table of Contents'. Questions concerning any part of the manual may be directed to [biosend@iu.edu](mailto:biosend@iu.edu) for further clarification.

[Manual of Procedures ↓](#)

### TRAINING SLIDES

These slides correspond to the BioSEND DxCTEII protocol training. Training is available upon request by contacting [biosend@iu.edu](mailto:biosend@iu.edu).

[Training Slides ↓](#)

### SAMPLE SHIPPING

BioSEND can receive samples Monday-Friday, excluding holidays. Frozen samples should be shipped M-W. Ambient samples may be shipped on Th.

[Generate UPS airbill or schedule pickup →](#)  
[Check holiday closures →](#)  
[What do I do for Friday blood draws →](#)

# Contacts

---

## Indiana University

General Questions/Shipment Notifications:

[biosend@iu.edu](mailto:biosend@iu.edu)

317-278-6158

Request kits:

<http://kits.iu.edu/biosend/dlb>