



Biospecimen Exchange for Neurological Disorders

Study in Parkinson Disease of Exercise Phase 3 Clinical Trial: SPARX3

BIOSPECIMEN COLLECTION & PROCESSING

Overview

1. Specimen uniformity and quality
2. Site Equipment
3. Procedures
 - Kit Contents and Ordering
 - Sample Labeling
 - Sample Collection and Processing
 - Shipping Samples
 - Non-Conformance
4. Contact Information

Specimen Uniformity and Quality

GENERAL REMINDERS

Specimen Standardization and Quality

Most biomarkers are sensitive to *time* and *temperature*

- Standardization of processing across sites is key
- Specimens must be processed within 2 hours of collection
- Reference the *BioSEND Biomarker Specimen Collection, Processing, and Shipment Manual* as needed
- Do not replace or supplement any kit components without first receiving approval from BioSEND/NINDS

Questions? Email biosend@iu.edu

Site Consumables and Equipment

Sites will need to supply the following items:

- Gloves
- Alcohol wipes
- Butterfly needles
- Tourniquet
- Gauze pads
- Bandages
- Microcentrifuge tube rack
- Sharps bin and lid
- Crushed ice
- Pipettes and pipette tips
- Centrifuge capable of maintaining 4°C
- -80°C Freezer
- Dry ice

Procedures

MAINTAINING SPECIMEN UNIFORMITY AND QUALITY

Biospecimen Collection Protocol

Visit & Sample Type	BL	6M	12M	18M	24M
Serum aliquots, 1.5ml	6	6	6	6	6
Plasma aliquots, 1.5 ml	6	6	6	6	6
Buffy Coat	2	2	2	2	2
Whole Blood, 3ml	2	2	2	2	2

Kit Contents and Ordering

- All sites will be sent a Supplemental Kit with their first kit shipment
 - Contains extra blood collection tubes and processing supplies
 - May be used to replace items in study visit kits
- Study Visit Kits should be ordered as soon as visits are planned
 - Contains collection, processing, and shipping supplies specific to each visit
 - Includes barcoded labels
 - The supplies/labels in each study visit kit are intended for that visit only


Kit Contents and Ordering – REDCap Survey

<http://kits.iu.edu/biosend/sparx3/>

Order kits online through the Kit Request Module for:

- Blood Kits
- Supplemental Kit
- Extra Supplies

Please provide as much notice as possible when ordering kits and/or supplies.



The screenshot shows the 'NINDS Biosend SPARX3 Kit Request System' interface. At the top, the 'BioSEND' logo is displayed in large blue letters, with the tagline 'Biospecimen Exchange for Neurological Disorders' below it. The title 'NINDS Biosend SPARX3 Kit Request System' is centered in blue. Below this, there is a 'Study Site' label and a dropdown menu. A red asterisk and the text '* must provide value' are positioned below the 'Study Site' label. At the bottom right, there is a 'Submit' button.

BioSEND
Biospecimen Exchange for Neurological Disorders
NINDS Biosend SPARX3 Kit Request System

Study Site

* must provide value

Submit

Kit Contents and Ordering: Confirm Site Info

SPARX3 Kit Request Module

Study Site <small>* must provide value</small>	Indiana University	Select your site from the drop-down list	
Indiana University School of Medicine Carolyn Dunifon Dept. of Medical & Molecular Genetics 351 West 10th Street, TK-318 Indianapolis, IN 46202-3002 (317) 274-5751 cdunifon@iu.edu			
Is the contact name above correct? <small>* must provide value</small>		<input type="radio"/> Yes <input checked="" type="radio"/> No	Verify contact information and update if needed
		reset	
New Contact Name <small>* must provide value</small>		Claire Wegel	
Is the shipping address above correct? <small>* must provide value</small>		<input type="radio"/> Yes <input type="radio"/> No	reset
Is the e-mail address above correct? <small>* must provide value</small>		<input type="radio"/> Yes <input type="radio"/> No	reset
Is the phone number above correct? <small>* must provide value</small>		<input type="radio"/> Yes <input type="radio"/> No	reset

Kit Contents and Ordering: Kit Types

SPARX3 Kit Request Module

- Kits are not specific to a subject or time point. After collection, sites will indicate the subject and time point to which BioSEND should link the samples.

Kit Type **Please allow two weeks for shipment** * must provide value	<div><input checked="" type="radio"/> Blood Collection Kit</div> <div><input type="radio"/> Supplemental Kit</div> <div><input type="radio"/> Extra Supplies</div> <p>Please specify in comments if you need kits before the standard two week shipment time.</p>
Blood Collection Kit Quantity * must provide value	<input type="text"/>
Comments	<div><input type="text"/></div> <div>Expand</div>

Kit Contents and Ordering: Kit Breakdown

SPARX3 Kit Request Module

Comments
<div><div></div><div>Expand</div></div>
<p>Each Blood Collection Kit contains:</p> <ul style="list-style-type: none">2 - Purple-top EDTA tube (glass, 10 ml)2 - Red-top Serum tube (glass, 10 ml)2 - Purple-top EDTA tube (plastic, 3ml)6 - Bubble-tube sleeve6 - Micronic cryovial (2ml) - purple6 - Micronic cryovial (2ml) - red2 - Micronic cryovial (2ml) - grey2 - Disposable transfer pipet (3ml)1- 48-slot micronic cryovial tray2 - Biohazard bag with absorbent sheet1 - Shipping label packet (Dry Ice, Fragile, UN3373)1 - Airway bill envelope1 - Dry Ice Shipper1- Specimen Label Set-- ST labels



Kit contents of selected kit will appear at the bottom of the page

Kit Contents and Ordering: Blood Kit



Kit and Supply Ordering

- Click “Submit” to send order to BioSEND; staff will confirm receipt of your order
- Please allow two week turn-around time for kit shipments
- If urgent request needed, please note date needed by in comments and email BioSEND. We cannot guarantee urgent orders, but we will do our best to accommodate.
- BioSEND will send confirmation of shipment and tracking number when supplies ship


Sample Labelling

Labels are provided by Indiana University for BioSEND samples only


- Please check that all samples are properly labelled with correct specimen type



Kit Labels



Specimen (Collection tube) Labels




Cryotube Labels

Sample Labelling


Kit Number

- Any kit number can be used for any subject or visit—which is why it is very important to document which kit number was used for which subject and visit


Each Blood Collection Kit will have a unique “Kit Number” used to track participant samples and provide quality assurance



Kit Labels



Specimen (Collection tube) Labels




Cryotube Labels

Sample Labelling – Kit Labels

Kit Number labels are placed:

- On the biohazard bag of the cryovial transport box and biohazard bag with EDTA.
- On the lid of frozen shippers



Sample Labelling – Collection Tube Labels



Study



Kit Number



Unique Barcode




Sample Type




Collection Tube

Sample Labelling – Collection Tube Labels



Sample Labelling – Cryotube Labels



Biorepository/Study



Kit Number




Sample Type

Sample Labelling – Cryotube Labels

Note: cryotube barcodes are printed on tube, not on the label. Barcodes are tied to the specific kit they arrive in (ie, a specific subject and visit). **Please do not mix the labels from one kit with the cryotubes from another.**

Barcode is in human readable format on side of tube, 2D format on bottom of tube



Sample Labelling – Cryotube Labels


Please verify that you are using the correct cryotube for the correct specimen type

Cap Color	Specimen Type
Purple	Plasma
Grey	Buffy Coat
Red	Serum

Sample Collection and Processing


Blood Tube Draw Order

2 x 10ml Serum




1

2 x 10ml EDTA
for plasma and
buffy coat




2

2 x 3ml EDTA
for whole
blood




3


Sample Collection and Processing: Serum



Sample Collection and Processing: Plasma & Buffy Coat





Sample Collection and Processing: Whole Blood



Sample Collection and Processing: Aliquots


Filling biomarker serum and plasma aliquots:

- Fill as many cryovials as possible to 1.5 ml (plasma & serum)
- Over-filled vials may burst in freezer!
- Ship ALL material to IU, even if final vial is less than standard volume



Sample Collection and Processing: Timeline

Timeline for blood processing



Sample Collection and Processing: Issue #1

Troubleshooting Blood Collection

Issue #1: Tube with little/no vacuum

- Always check expiration date on the tube before beginning blood draw and discard expired tubes
- Store tubes at “room temperature” – extreme temperature can affect vacuum
- Keep extra vacutainer tubes from supplemental kit nearby during blood draw to replace “bad” tubes
- If this is a frequent occurrence, report tube type and lot number to IU.

Sample Collection and Processing: Issue #2

Troubleshooting Blood Collection

Issue #2: Hemolyzed serum and/or plasma caused by incorrect collection

Cause: Blood Collection Methods	Corrective Action
Improper venipuncture site	Draw from median cubital, basilic, and cephalic veins from antecubital region of arm
Prolonged tourniquet use	Tourniquet should be released after no more than 1 min, excessive fist clenching should be avoided
Not allowing alcohol to dry on skin before venipuncture	Without touching, allow the venipuncture site to air dry
Use of too large/small bore needle resulting in excess force applied to blood	Avoid using too small/large needle. Needle size dependent on the subject's physical characteristics & amount of blood to be drawn. Most commonly used sizes are 19 – 23.
Pulling/pushing plunger too fast while drawing/transferring blood	Avoid drawing the syringe plunger too forcefully when collecting blood
Ensure all blood collection assemblies are fitted securely, to avoid frothing	

For more information, visit: http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk_Jan2004_VS7167.pdf

Sample Collection and Processing: Issue #2 continued

Troubleshooting Blood Collection

Issue #2: Hemolyzed serum and/or plasma caused by incorrect processing

Cause: Sample Processing Methods	Corrective Actions
Vigorous mixing/shaking	Gently invert blood collection tube when mixing additive with specimen, follow guidelines in Biologics Manual regarding number of times to invert each type of tube
Not allowing serum to clot for recommended time	Serum tubes without clot activator should be allowed to clot for 60 min in a vertical position
Exposure to excessive heat or cold	Keep samples at ambient temperature until processing
Prolonged contact of serum/plasma with cells	Do not store uncentrifuged samples beyond recommended time

For more information, visit: http://www.bd.com/vacutainer/pdfs/techtalk/TechTalk_Jan2004_VS7167.pdf

Sample Collection and Processing Form

Direct link:

<https://redcap.link/SPARX3SampleForm>

First page captures basic subject and visit information



Biospecimen Exchange for Neurological Disorders

Please complete the Specimen Collection and Processing Form, below.

 Returning?

A A A



Page 1 of 2

BioSEND SPARX3 Study

Study Site

Email address of staff member completing this form

Note: A copy of the completed sample form and the shipping manifest will be sent to this address.

GUID:

Sex (used for DNA quality control)

Visit

IU Kit Number

6 characters remaining

Next Page >>

Save & Return Later

Sample Collection and Processing Form

Direct link:

<https://redcap.link/SPARX3SampleForm>

Second page captures processing information

Blood Collection and Processing

Date of venipuncture blood collection

 Today M-D-Y

Time of venipuncture blood collection

 Now H:M

Use 24 Hour clock

Date participant last ate

 Today M-D-Y

Time participant last ate

 Now H:M

SERUM (red-top tubes, 10 mL)

Was blood collected and processed for SERUM?

Yes

No

reset

Time of SERUM tube centrifugation

 Now H:M

Use 24 Hour clock

Duration of SERUM tube centrifugation

15

minutes

Rate of SERUM tube centrifugation

1500

x g

Temperature of SERUM tube centrifugation

4

degrees Celsius

Total volume of SERUM collected

ml

Sample Collection and Processing Form

Direct link:

<https://redcap.link/SPARX3SampleForm>

PDF form of responses will be emailed to you. **Print a copy of the Frozen Shipping Manifest and include with shipment.**

SPARX3 Frozen Shipping Manifest

Page 1

Please verify/update the information below. When you click the "Submit" button below, a PDF copy of the Frozen Shipping Manifest will be emailed to you for Subject [subj_id].

Please print a copy of that document and include it in the Kit #[kit_num] shipping container.

Study Site:

- ☐ Boston University (Charles River Campus)
- ☐ Cleveland Clinic
- ☐ Columbia University Medical Center
- ☐ Emory University
- ☐ Iowa State University
- ☐ Kent State University
- ☐ Louisiana State University
- ☐ Mayo Clinic, Rochester
- ☐ New York University Langone Health
- ☐ Northwestern University
- ☐ Ohio Health
- ☐ Oregon Health & Science University
- ☐ Rush University Medical Center
- ☐ University of Alabama at Birmingham
- ☐ University of Alberta
- ☐ University of California, San Francisco
- ☐ University of Cincinnati
- ☐ University of Colorado, Denver
- ☐ University of Florida
- ☐ University of Michigan
- ☐ University of Minnesota
- ☐ University of Pennsylvania
- ☐ University of Pittsburgh
- ☐ University of Southern California
- ☐ University of Texas Medical Branch
- ☐ University of Texas San Antonio
- ☐ University of Utah
- ☐ University of Virginia
- ☐ Washington University St. Louis
- ☐ Wildred Laurier University of Canada

GUID:

Visit:

- ☐ BL
- ☐ 6M
- ☐ 12M
- ☐ 18M
- ☐ 24M

IU Kit Number:

Date of blood collection:

Sample Collection and Processing Form

Submission of the Shipping Manifest portion of the form serves as shipment notification to BioSEND

- Must be completed prior to shipment
- If samples are shipped and not received, BioSEND will follow-up with courier. It is recommended that sites also track shipment to ensure safe delivery

Shipping Frozen Samples: Tips

Packing and Shipping Frozen Samples

- SPARX3 samples all ship frozen on dry ice
- Frozen samples should be shipped **only** Monday through Wednesday
- Always fill carton to **top** with dry ice
- Do not pack shipment until the day of pickup
- Do not overlap the dry ice label with any other labels or the airway bill




Shipping Samples: Frozen

Do not ship more than 4 biohazard bags in a single shipper (equivalent to two subject-visits).

- Allows room for dry ice to keep samples frozen in transit
- Minimizes loss in the event of courier issues



Shipping Samples – UPS: <https://kits.iu.edu/UPS>



INDIANA UNIVERSITY

SCHOOL OF MEDICINE

and Affiliated Biorepositories

Search for address

Code

Company

Contact

Address 1

Address 2

Address 3

City

State/Province

Postal Code

Country/Territory

Icahn School of Medicine - Mt. Sinai

Kenny Persaud

1425 Madison (Icahn Building)

Icahn L4 39

New York

NY

10029

United States

Shipment Information

Study Group

Weight

Dry Ice Weight

Description of Return

SSBC

20

10

Biologic Specimens

Choose Study

Enter weight

Click "Ship"

Reset

Ship

Shipping Samples via UPS

IU UPS ShipExec Shipping Portal

- Print out UPS air waybill
- Ensure all elements (barcode, return address, etc.) printed clearly
- Fold and insert UPS air waybill into clear plastic sleeve on package

JOHN SMITH INDIANA UNIVERSITY 410 WEST 10TH STREET INDIANAPOLIS IN 46202	2 LBS	1 OF 1
RS		
SHIP TO: SCHOOL OF MEDICINE 317-278-2694 INDIANA UNIVERSITY TK 217 351 W 10TH ST INDIANAPOLIS IN 46202		
	IN 461 9-01	
		
UPS NEXT DAY AIR		1
TRACKING #: 1Z 976 R8W 84 3985 8595		
 SAMPLE		
BILLING: P/P DESC: Biological Specimens RETURN SERVICE		
Reference No.1: 4087277		
XOL 20.03.09 NV45 83.0A 12/2019		 TM

Shipping Samples: Closures

Date	Holiday
January 1	New Year's Day
3 rd Monday in January	Martin Luther King, Jr Day
4 th Monday in May	Memorial Day
June 19	Juneteenth (observed)
July 4	Independence Day (observed)
1 st Monday in September	Labor Day
4 th Thursday in November	Thanksgiving
4 th Friday in November	Friday after Thanksgiving
December 25	Christmas

Please also consider weather when shipping. UPS will post service updates on their webpage. Reach out to BioSEND if you are unsure if it is safe to ship.

Non-Conformance

Non-conformance to standard procedures may reduce the utility of the biospecimens:

- Not processing serum/plasma within 2 hours of collection allows for breakdown of certain proteins and small molecules
- Over/under centrifuging changes plasma, serum composition



Non-Conformance Reporting con't

Most common non-conformance issues:

- Shipment notification not sent
- Samples shipped for weekend/holiday delivery
- Sample form incomplete/inaccurate
- Low volume
- Unlabeled or mislabeled tube(s)
- Sample hemolysis



Non-Conformance and Inventory Reporting

Most common non-conformance issues:

- BioSEND will notify sites directly of any issues upon receipt
- BioSEND will email sites a monthly inventory report of all samples received from that site to date
- ***If you are experiencing issues, please reach out to us for help! It is much easier to prevent an issue before sample collection & shipment than trying to fix it after the fact***

Contacts

Indiana University

General Questions/Shipment Notifications:

biosend@iu.edu

317-278-6158

Request kits:

<http://kits.iu.edu/biosend/sparx3/>